


FIG.4









Inventor: CHU et al.

Docket No.: 7432.141US01
Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY
Serial No.: 09/804866

Sheet 7 of 17







Serial No.: 09/804866

## FIG.7A

## **Option Explicit**

``` ///. Type Declaration ``` **////.** Constant Declaration Const ciFrmHeight As Integer = 5900 ////. Form height Const ciFrmHeight As Integer = 9600 ///. Form height **//// Variable Declaration** Dim icountFound As Integer ///. No of series found ///. Prototypes ```` /// Funtion Declaration Private Sub Form Load() ///. Initialization Call P LoadForm End Sub



Inventor: CHU et al. Docket No.: 7432.141US01

Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY

Serial No.: 09/804866

Sheet 9 of 17

## FIG.7B

```
```
////.
///. Synopssis
             :P_LoadFoam()
///, Name
             :Private Sub P_LoadForm()
////.
///. Input
             :None
////.
///. Return
             :None
Private Sub P_LoadFoam()
 ///, Label initialization
 IblFirstRow.Caption = " "
 lblFirstValue.Caption = " "
 IblSecondRow.Caption = " "
 IblSecondValue.Caption = " "
 IbIThirdRow.Caption = " "
 IbIThirdValue.Caption = " "
 /// Parameters initialization
 'glErrCode = RC_SUCCESS
 'gsErrMsg = " "
End Sub
```



Serial No.: 09/804866 Sheet 10 of 17

Private Sub cmdOK\_Click()
Call P\_FirstUpdate

```
End Sub
```
////.
////. Synopssis :P_FirstUpdate
////.
///. Name
             :Private Sub P_FirstUpdate()
////.
//// Description: Update in the first form
////.
////. Input
             :None
////.
///. Return
             :None
////.
Private Sub P FirstUpdate()
 Screen.MousePointer = vbHourglass
 'glErrCode = RC_SUCCESS
 Screen.MousePointer = vbDefault
 Call FindNoSeries
End Sub
Private Sub cmdCancel_Click()
 Unload Me
 Exit Sub
End Sub
Private Sub Form_Unload(Cancel AS Integer)
 ///. Clean up the form
 Call P UnloadForm
End Sub
```



Docket No.: 7432.141US01 Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY

Serial No.: 09/804866

```
////.
////. Synopssis :P FirstUploadForm
////.
///. Name
            :Private Sub P UploadForm()
////.
//// Description :Release the resource which occupied by the parameters
////.
            in the form
////, Input
             :None
////.
///. Return
            :None
////.
```
Private Sub P_UploadForm()
End Sub
```
////.
////. Synopssis :FindNoSeries
////.
///. Name
            :FindNoSeries() As Long
////.
//// Description :Function to update Uex file and stock cosing price
////.
///. Input
             :N/A
////.
///. Return
            :FindNoSeries RC SUCCESS
  Success
////.
                  RC_ERROR Error
////.
Public Funtion FindNo Series() As Long
 Dim gsPath As String
 Dim gsCEFile As String
```

Inventor: CHU et al. Docket No.: 7432.141US01

Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY

Serial No.: 09/804866

Sheet 12 of 17

Dim gsErrMsg As String

Dim gsExist As String

Dim takeFile As String

Dim giCEFileNo As Integer

Dim bCloseFileNo As Boolean

Dim sNextLn As String

Dim asNumberFound(800) As String

Dim sFileName As String

'FindNoSeries = RC\_ERROR

bCloseFile = True

'gsErrMsg = " "

'gsPath = App.Path

////, Name of file is SortRecord

gsCEFILE = gsPath & "/" & "SortRecord" & "/"

sExit = Dir(gsCEFile)

sFileName = txtFileName.Text & ".txt"

If sExist = sFileName Then

bCloseFile = False

End If

If bCloseFileThen

////. "not found. Update process not proceeded."

gsErrMsg = "SortRecord not found. Process stops."

MsgBox (gsErrMsg)

'FindNoSeries = RC\_ERROR

**Exit Function** 

End If

giCEFileNo = FreeFile

takeFile = gsCEFile & sExist

Open takeFile For Input As #giCEFileNo

Dim tlLnNum As Long



Docket No.: 7432.141US01

Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY

```
ILnNum = 0
Line Input #giCEFileNo. sNextLn
Do Until UCase(Trim(sNextLn)) = "END"
 ///. Vailidation
 If EOF(giCEFileNo) Then
   ///, "Unexpected end of file."
   gsErrMsg = "Unexpected end of file."
   MsgBox (gsErrMsg)
   'FindNo Series = RC ERROR
   Exit Function
 End If
 If NotIsNumeric(sNexrLn) Then
   ///, "Unexpected end of file."
   gsErrMsg = "Non-numeric input" & "Row" & ILnNum + 1 & "."
   MsgBox (gsErrMsg)
   'FindNoSeries = RC_ERROR
   Exit Function
 End If
 If ILnNum <> 0 Then
  If Val(asNumberFound(ILnNum - 1)) < Val(sNextLn) Then
    gsErrMsg = "The input numbers are not sorted.(at row" & ILnNum + 1 & ".)"
    MsgBox (gsErrMsg)
   'FindNoSeries = RC ERROR
    Exit Function
  End If
```

**DoEvents** 

End If

Line Input #giCEFileNo. SnextLn

asNuberFound(ILnNum) = SNextLn

ILnNum = (ILnNum + 1)

Loop



Sheet 14 of 17

## If Not EOF(giFileNo) Then

///. "Unexpected end of file"
gsErrMsg = "unexpected end of file"
MsgBox (gsErrMsg)
'FindNoSeries = RC\_ERROR
Exit Function

End If

Close #giCEFileNo

Dim IPrimaryNo As String Dim LsecondaryNo As String Dim ITotalNo As String Dim iOuterLoop As Integer Dim ilnnerLoop As Integer Dim bFound As Boolean Dim bNotFound As Boolean Dim IMagicNo As String Dim ITolNo As String Dim IMagicOne As String Dim IMagicTwo As String Dim ITolOne As String Dim ITolTwo As String Dim asNumberOuput(200) As Long Dim iCountOutput As Integer Dim iFileNum As Integer Dim putFile As String

bCloseFile = True

///. "Unexpected end of file"

'gsPath = App.Path

///. Name of file is SortRecord

gsCEFile = gsPath & "/" & "Result" & "/"



Docket No.: 7432.141US01 Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETRY Serial No.: 09/804866

```
'sExist = Dir(gsCEFile)
"If sExist = "Try I .txt"
  bCloseFile = False
'End If
If bCloseFileThen
 //// "not found. Update process not proceeding."
 gsErrMsg = "Result not found. Process stops."
  MsgBox (gsErrMsg)
  'FindNoSeries = RC ERROR
  Exit Function
End If
///. Open file for writing result
iFileNum = FreeFile
putFile = App.Path & "\" & "Result & "\" & txtFileName.Text & ".txt"
Open putFile For Append as #iFileNum
///, Set initial value
iCountFound = 0
clMagicOne = Vai(txtDiff1.Text)
clMagicTwo = Vai(txtDiff2.Text)
clTolOne = Vai(txtTol1.Text)
clTolTwo = Vai(txtTol2.Text)
bFound = False
bNotFound = False
iCountOuput = 0
iMagicNo = clMagicOne
iTotalNo = ILnNum
For iOuterLoop = 0 To (ITotalNo - 1)
  IPrimaryNo = asNumberFound(iOuterLoop)
  ///, Set all captions
  lblFirstValue.Caption = IPrimaryNo
  IblFirstRow.Caption = IiOuterLoop = 1
```



Title: SEQUENCING OF PEPTIDES BY MASS SPECTROMETR Serial No.: 09/804866 Sheet 16 of 17

```
IblSecondValue.Caption =
IblSecondRow.Caption = " "
IblThirdValue.Caption = " "
IbIThirdRow.Caption = " "
///. Reset variable in interloop
   bFound = False
   bNotFound = False
   iMagicNo = clMagicOne
   iTolNo = ILnNum
////. set interloop begin value
ilnnerloop = iOuterLoop = I
///, do until end of list or pattern found or not found
Do until BFound Or bNotFound Or (ilnnerloop > ITotal - 1)
   ISecondaryNo = asNumbersFound(ilnnerLoop)
   If (Val(PrimaryNo) - Val(SecondaryNo)) .= (Val(MagicNo) - Val(ITolNo)) And_
   (Val(PrimaryNo) - Val(SecondaryNo)) .= (Val(MagicNo) - Val(ITolNo)) And_
      If IMagicNo = clMagicTwo Then
      IbIThirdValue.Caption = ISecondaryNo
      IbIThirdRow.Caption = iInnerLoop + 1
      bFound = True
      'FindNoSeries = RC SUCCESS
      Write #iFileNum, IPRimaryNo & "&" & ISecond & _
          "&" & ISecondaryNo
      Write #iFileNum, IPRimaryNo & "&" & ISecond & "&" & ISecondaryNo
         '((IPrimaryNo - 18 + ISecond + ISeconaryNo +28) / 3
      iCountFound = iCountFound +1
   Else
      lblSecondValue.Coupon = ISecondaryNo
      lblSecondRow.Coupon = ilnnerLoop + 1
      ISecond = ISecondaryNo
      IMagicNo = clMagicTwo
    End If
Elself (Val(SecondaryNo) - Val(PrimaryNo)) > Val(MagicNO) Then
```



BNotFound = True

End If

iInnerLoop = iInnerLoop + 1

Next

///. Close file

Close #iFileNum

If bFound = False Then

'FindNoSeries = RC\_NOT\_FOUND

End If

If iCountFound <> 0 Then

IbITrueFalse.Caption = "True"

MsgBox (iCountFound & "specific series number are found and save in

"&" txtFileName.Text & ".txt)

Else

IblTrueFalse.Caption = "False"

MsgBox ("The specific series number is not found!!!")

End If

**End Function**